

magstim[®] TMS THERAPY

HORIZON[®]

Performance



Help your patients **REDISCOVER THE COLOUR OF LIFE**

DELIVERING PERFORMANCE, VERSATILITY AND RELIABILITY

As The Brains Behind TMS™, we believe in delivering high performance, versatile, TMS solutions so that you can provide the best treatment for your patients.


The HORIZON® performance package combines a comprehensive, all-in-one TMS solution with our 30 years of expertise to provide a versatile, complete TMS package.

HORIZON® Performance

- ✓ **Performance** | The HORIZON® energy recovery system delivers consistent stimulation throughout treatment.
- ✓ **Effortless** | The E-z Arm provides effortless manipulation of coil placement and secure coil location.
- ✓ **Versatile** | Theta burst & standard TMS capabilities, with the option to save personalised patient protocols.
- ✓ **Reliable** | Developed by the pioneers who designed TMS. With over 30 years experience, you can rely on our expertise.
- ✓ **Convenient** | A comprehensive all-in-one solution, the E-z arm provides for weightless, secure and versatile coil positioning.


No Pulse Decay

Consistent stimulation during therapy



Versatility:

inc. 19 minute protocols & Theta Burst (TBS)



WHAT DO THE ROYAL COLLEGE OF PSYCHIATRY SAY ABOUT TMS?

“Over the last decade rTMS has been used widely for the treatment of depression and is now an established safe and effective treatment option for depression and treatment resistant depression.” Source: Royal College of Psychiatrists statement on rTMS for depression. Position statement CERT03/17



E-z Coil Arm

Secure, safe and versatile, the E-z arm provides full maneuverability of the coil which can be held securely at almost any position.

E-z Cool Coil

The intelligent cooling system automatically adjusts to the temperature of the coil providing uninterrupted protocols and a faster patient throughput.

HORIZON® Stimulator

The HORIZON® energy recovery system ensures consistent performance in every session.

High powered and theta burst capabilities, easy plug-and-go use from the pre-programmed protocol and the versatility of saving custom patient protocols.

HORIZON® MT Remote Coil

Ergonomic, hand-held coil with integrated control to simplify motor threshold assessment



Therapy Chair

Mobile and motorized, for ease of use, the therapy chair is fully adjustable and secure for optimal patient and operator comfort.



E-z Cart

Convenient, highly mobile and secure. The E-z cart provides an all-in-one solution making portability, safety and storage simple.



magstim® TMS THERAPY


For more information, or to request a quote; call: **+44 (0) 1994 240798**

www.magstim.com | salesoperations@magstim.com




The Magstim Company Ltd.
Whitland Industrial Estate, Spring Gardens
Whitland, Carmarthenshire
SA34 0HR, UK

MK1806-02



Depressed Brain



After TMS

What is TMS?

The majority of depression episodes (>63%) result in non-remission after first-line treatment with any chosen anti-depressant medication¹.

Transcranial Magnetic Stimulation (TMS) involves a series of repetitive, brief and highly focused magnetic pulses used to stimulate brain cells. TMS induces a small electric current in the targeted area of the brain, increasing brain activity to normal levels.

HORIZON™ TMS Therapy is non-invasive, non-systemic and has few known side effects²; resulting in better patient outcomes.

1. Rush A J et al. Acute and longer-term outcomes in Depressed outpatients requiring one or several treatment steps: A STAR*D Report. Am J Psychiatry. 2006. 163. 1905-1917
2. Carpenter, L et al. Transcranial Magnetic Stimulation (TMS) for Major Depression: A multi-site, naturalistic, observational study of acute treatment outcomes in clinical practice (2012) J Depress Anxiety 2012; 29, 587-596.

CONTACT US TODAY AND START CHANGING LIVES