

Feature	Horizon® 3.0	Horizon® 3.0 StimGuide® +	Horizon® 3.0 StimGuide® + Magstim® Connect
Standard rTMS and iTBS protocols	✓	✓	✓
Option for tailoring protocols	✓	✓	✓
Proprietary algorithm delivers consistent stimulation throughout treatment with no pulse decay	✓	✓	✓
Ergonomic and secure coil positioning with enhanced braking system – arm stabilization	✓	✓	✓
Increased coil lifetime	✓	✓	✓
Modern, intuitive, adjustable touchscreen user interface (UI)	✓	✓	✓
Efficient smart cooling capabilities for maximum patient output	✓	✓	✓
Trigger in/out capabilities	✓	✓	✓
Multi-layer real-time quality monitoring		✓	✓
Ez-MT®, 2-step motor threshold determination with EMG option		✓	✓
FDA and UKCA approved 3D navigation system for accurate, precise and consistent coil placement		✓	✓
Real-time XYZ tracking with four parameters – distance, contact, tilt and rotation		✓	✓
Powerful quality and pulse quality scores and treatment analytics		✓	✓
FDA registered medical device data system*			✓
Downloadable patient response reporting for insurance and referral providers			✓
Remote assignment of protocol and treatment templates			✓
Mobile depression surveys and score input			✓
Synchronized data between multiple devices to support multi-site clinics.			✓


Connected Care from Anywhere.



Magstim®
 The original TMS architects continue to innovate, bringing you the Horizon 3.0 with StimGuide+ and Magstim Connect, to further enhance your TMS program and patient outcomes.

magstim HORIZON® 3.0

Quality Control. Consistency. Confidence.



Imagine the most advanced clinical TMS system that puts you in control — to deliver customized, repeatable treatments, every time.

magstim
HORIZON® 3.0

magstim
STIMGUIDE® +

magstim
CONNECT

MAGSTIM.
A WELCONY BRAND

www.magstim.com +44 (0) 1994 240798

MK2101-02-UK

*Magstim Connect is a Medical Device Data System that is exempt from the pre-market notification procedures as described in 21 CFR § 807, Subpart E. Magstim Connect is listed with the US FDA.



H3 is a new class of TMS technology that enhances clinical workflows with navigated TMS and added layers of analytics — in an intuitive and versatile system

Meet the Horizon 3.0

- Pulse quality ‘score’ dashboard
- Mobile accessible clinical rating scales
- Patient specific protocols
- Referral and patient response reporting
- Connected care at all practice locations

Magstim Legacy


Magstim developed the first TMS system over 30 years ago and close to 16,000 studies published worldwide cite the use of Magstim stimulators. Magstim holds seven 510(k) clearances for its TMS therapy systems including StimGuide®, the first FDA-cleared 3D navigation system designed for the clinical setting. With the Horizon 3.0 with StimGuide+ clinicians can be confident they are delivering the right dose to the right location, every treatment, every time.

“(The H3) is fantastic. It really is so simple the way it is set up. It’s very patient friendly (and) clinically effective.”



— John Cavanagh, RestoreBrain

Delivering the right dose with every pulse

Magstim’s Proprietary Algorithm ensures you deliver correct and consistent dosage throughout a treatment session — without pulse decay — to optimize outcomes.



magstim STIMGUIDE+ Pulse Quality Score



magstim: CONNECT

Ensuring consistency across multiple networked locations

- Powerful patient data management and treatment reporting
- Connected care with networked functionality to support practice flexibility and remote upgrades
- Synchronized data between multiple devices in the HIPAA compliant cloud to support multi-site clinics
- Mobile depression surveys and score input

Magstim TMS Therapy is FDA and UKCA cleared for the indications of Major Depression Disorder, Anxious Depression, also known as Depression with Anxiety Symptoms and Obsessive Compulsive Disorder in adult patients. Magstim TMS is an effective, non-invasive, outpatient treatment, with few known adverse effects.

*Carpenter, L et al. Transcranial Magnetic Stimulation (TMS) for Major Depression: A multi-site, naturalistic, observational study of acute treatment outcomes in clinical practice (2012)” *Depress Anxiety* 2012; 29, 587-596.

Quality Control. Consistency. Confidence.